



***Insight* — Technical Overview 2.03**

Guidelines for Dielectric Cure Monitoring

Introduction

Following are some guidelines for successful dielectric cure monitoring. Although the method of DEA is simple—place a sample on the sensor, heat and start data acquisition—care and good practice are very helpful for making good measurements.

Preparing the sensor

- Do remove oils or other contaminants by cleaning sensors with acetone, alcohol or other solvent
 - Remove adsorbed solvents, which might interfere with test measurements in air, by heating sensors above 100 °C for a short time
- Do prevent material from adhering to reusable sensors by applying silicone or non-conductive mold release, as shown in Figure 1





Figure 1
Apply mold release to reusable sensors

- Do prevent short circuits by avoiding contact between sensor bond pads and electrically conductive surfaces
 - Cover bond pads with Kapton® or polyimide tape, as shown in Figure 2



To prevent short circuit of bond pads
—Avoid electrically conductive surfaces in this area
—Insulate above and below bond pad area

Figure 2
Insulating bond pad area

- Do run sensor leads parallel to each other to reduce capacitance between leads
 - Don't use twisted leads, which increase cable or base capacitance

Laying up the sample

- Do place sample over entire electrode area, as shown in Figure 3




Figure 3
Sample application area on sensor

- Do use samples thicker than the separation between interdigitated electrodes
 - When the sample is thinner than the electrode separation, the sensor will also detect air or material on the top side of the sample

- Do stack at least two or three layers of prepreg on top of a sensor, to ensure enough resin for good measurements
- Do use a filter between the sensor and materials containing graphite or conductive fibers, as shown in Figure 4
 - Glass cloth with small pore size, fiberglass felt or laboratory filter paper are good filter materials
- Do use aluminum foil or release film above and below the lay-up, as shown in Figure 4, to prevent the sample from adhering to platen or mold surfaces
- Do apply pressure to solid samples, or solid samples that melt during processing, for good contact with the sensor




Figure 4
Suggested lay-up for prepregs

Reducing measurement noise

- Don't use long, unshielded leads, which act like antennas and can pick up interference
 - Signal levels at the end of cure are low and measurements are more susceptible to noise at this time
 - Use coaxial cable with guarded or grounded shields
- Don't place the sensor on or near large ungrounded metal surfaces
 - Ungrounded metal acts like an antenna that can pick up noise, which the sensor detects
 - Ground metal surfaces around the sensor
- Don't place power cords near the sensor
 - AC mains voltages are around 100 – 240 VAC but sensor signals may be in the range of only 10 mV




Figure 5
Ways to reduce noise in dielectric sensors

Interpreting dielectric data

- Do look at loss factor, as in Figure 6, when studying the dipole response
 - Total loss factor is the sum of loss from the flow of mobile ions and the loss from dipole rotation
 - Using loss factor makes it easier to see and separate the mobile ion response from the dipole response
 - Loss factor is inversely proportional to frequency when mobile ions dominate response




Figure 6
Loss factor of curing polyurethane

- Do look at resistivity from multiple frequencies, as in Figure 7, to determine ion viscosity
 - Ion viscosity indicates cure state
 - Proportional to change in viscosity before gelation
 - Proportional to change in modulus after gelation
 - Ion viscosity is frequency independent resistivity due to mobile ions
 - Ion viscosity dominates where curves from different frequencies closely or completely overlap
 - Dipole response dominates where curves from different frequencies do not overlap




Figure 7
Resistivity of curing polyurethane (plotted against ion viscosity axis)

- Do use ion viscosity from a single, optimum frequency—if possible—to calculate slope, as shown in Figure 8
 - Lower frequencies show ion viscosity response better for whole cure
 - Higher frequencies show dipole response better at end of cure




Figure 8
Ion viscosity and slope for curing polyurethane

Conclusion

With proper sample preparation, lay-up, shielding of leads and attention to the electrical environment it is possible to make good, reproducible measurements of dielectric properties for monitoring the cure of thermoset and composite materials.



Lambient Technologies LLC
649 Massachusetts Avenue, Cambridge, MA 02139
857-242-3963
www.lambient.com
info@lambient.com